User manual

PrimeFORCE

G5 series

Language index

English: Page 3-37

Nederlands: Pagina 38-72

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Français: Page 108-142

Table of contents

Introduction	4
Red light therapy	5
Specifications	6-10
Benefits of Liroma red light therapy	11
Instructions for use	12
Product structure	13
What's included?	14
Accessories	15
Installing the door hanger	16-17
MULTI setup: connecting two Elite RLT panels	18-19
Touchscreen control panel	20
Settings	21-25
Red light therapy presets	26
How does it work?	27
Understanding wavelengths	28-30
RLT Treatment Protocols	31-33
Warnings	34-36
Care and cleaning	36
About Liroma	37

Thank you for your purchase!

Thank you for choosing Liroma red light therapy!

Liroma wants its customers to feel healthy, fit and relaxed. To achieve this, Liroma has a wide range of high quality red light therapy products. Our products help with complaints such as: fatigue, poor circulation, pain in muscles and joints, pain due to underlying diseases, poor sleep, skin problems, stress and many other ailments. Our products target both mental and physical ailments.

Note: Red light therapy lamps do not get hot! During treatment, however, you may feel heat and sweat inside.

Note: Near-infrared light is invisible to the naked eye. It may appear as if some LEDs are not working, but this is completely normal.

Your wellness journey starts here!

Settle in, enjoy the moment, and let the soothing glow of red light support your body and mind. Here's to feeling stronger, sleeping better, and living each day with more vitality.

With love,

Red light therapy

Red light therapy uses a combination of red light and near-infrared light (NIR) to stimulate and repair the body at a cellular level. This therapy works by allowing the light to penetrate deep into the skin, with red light targeting the upper skin layers and NIR light reaching deeper into muscles and joints. A key effect of this light therapy is the stimulation of mitochondria—the energy factories of our cells—helping them produce energy more efficiently in the form of ATP (adenosine triphosphate).

With this increase in cellular energy, cells can recover more quickly, leading to reduced inflammation and improved blood circulation. Red light also promotes collagen production, which can contribute to firmer, more youthful skin and improved skin texture. Thanks to the combination of red light and NIR light, users often experience benefits such as reduced muscle and joint pain, enhanced skin health, faster muscle recovery, and increased energy levels. These effects make red light therapy a popular choice for those looking to support their well-being naturally.

SCAN and learn how red light therapy works

Specifications


PrimeFORCE Compact

Model	PrimeFORCE Compact / G5-100
Number of LEDs	72 5W Dual chip LEDs
Wavelengths (nanometer)	630nm · 650nm · 660nm · 670nm · 810nm · 830nm · 850nm · 1060nm
Light spectrum	Red light / Near-infrared light (NIR)
Power consumption	120 W
Input voltage	100V-240V 50/60Hz
Pulse frequency	1-10,000 Hz
Irradiance (15 cm)	163 mW/cm ² solar meter, 82 mW/cm ² spectrometer
Size	32 × 23 × 6.5 cm
Weight	3.9 kg
Beam angle	30°
Cooling fans	1
LED lifespan	100.000 hours expected LED lifetime
EMF emission	0,0 UT
Warranty	3 years
Best for	Targeted treatment

Specifications


PrimeFORCE Advanced

Model	PrimeFORCE Advanced / G5-200
Number of LEDs	144 5W Dual chip LEDS
Wavelengths (nanometer)	630nm · 650nm · 660nm · 670nm · 810nm · 830nm · 850nm · 1060nm
Light spectrum	Red light / Near-infrared light (NIR)
Power consumption	240 W
Input voltage	100V-240V I 50/60Hz
Pulse frequency	1-10,000 Hz
Irradiance (15 cm)	164 mW/cm ² solar meter, 88 mW/cm ² spectrometer
Size	64 × 23 × 6.5 cm
Weight	6 kg
Beam angle	30°
Cooling fans	3
LED lifespan	100.000 hours expected LED lifetime
EMF emission	0,0 UT
Warranty	3 years
Best for	Half body treatment

Specifications

PrimeFORCE Elite


Model	PrimeFORCE Elite / G5-300
Number of LEDs	288 5W Dual chip LEDs
Wavelengths (nanometer)	630nm · 650nm · 660nm · 670nm · 810nm · 830nm · 850nm · 1060nm
Light spectrum	Red light / Near-infrared light (NIR)
Power consumption	480 W
Input voltage	100V-240V 50/60Hz
Pulse frequency	1-10,000 Hz
Irradiance (15 cm)	162 mW/cm ² solar meter, 90 mW/cm ² spectrometer
Size	92 × 30 × 6.5 cm
Weight	11.1 kg
Beam angle	30°
Cooling fans	5
LED lifespan	100.000 hours expected LED lifetime
EMF emission	0,0 UT
Warranty	3 years
Best for	2/3 body treatment

Specifications

PrimeFORCE Multi


2 x PrimeForce Elite

Model	PrimeFORCE MULTI / G5-500
Number of LEDs	576 5W Dual chip LEDs
Wavelengths (nanometer)	630nm · 650nm · 660nm · 670nm · 810nm · 830nm · 850nm · 1060nm
Light spectrum	Red light / Near-infrared light (NIR)
Input voltage	100V-240V 50/60Hz
Pulse frequency	1-10,000 Hz
Irradiance (15 cm)	162 mW/cm ² solar meter, 90 mW/cm ² spectrometer
Size	184 × 30 × 6.5 cm
Weight	22.2 kg
Beam angle	30°
Cooling fans	10
LED lifespan	100.000 hours expected LED lifetime
EMF emission	0,0 UT
Warranty	3 years
Best for	Full body treatment

Specifications

PrimeFORCE Ultra

Model	PrimeFORCE ULTRA / G5-1000
Number of LEDs	1152 5W Dual chip LEDs
Wavelengths (nanometer)	630nm · 650nm · 660nm · 670nm · 810nm · 830nm · 850nm · 1060nm
Light spectrum	Red light / Near-infrared light (NIR)
Power consumption	1549 W
Input voltage	100V-240V I 50/60Hz
Pulse frequency	1-10,000 Hz
Irradiance (15 cm)	164 mW/cm ² solar meter, 85 mW/cm ² spectrometer
Size	189 × 58 × 6.4 cm
Weight	44.2 kg
Beam angle	30°
Cooling fans	36
LED lifespan	100.000 hours expected LED lifetime
EMF emission	0,0 UT
Warranty	3 years
Best for	Full-body treatment

Benefits of Liroma red light therapy

Product benefits

- Full-body or targeted treatment options
- Eight scientifically researched wavelengths
- Combination of red and near-infrared light
- Adjustable brightness levels
- Adjustable pulse frequencies (1–10,000 Hz)
- Touchscreen and remote control operation
- High irradiance for efficient treatments
- Flicker-free, non-invasive technology
- Suitable for home, wellness, and professional use
- No consumables or replacement parts required

Health benefits

- Supports healthy skin and collagen production
- Helps improve skin texture and appearance
- Supports muscle recovery after exercise
- Helps reduce post-workout muscle soreness
- Supports healthy circulation
- Supports joint comfort and mobility
- Helps promote cellular energy production (ATP)
- Supports the body's natural recovery processes
- May help reduce temporary inflammation
- Supports overall wellness and vitality

Instructions for use

Power switch: Turn the device ON/OFF via switch

Ideal usage frequency: Use the device 3–7 times per week. Consistent use generally provides the best results.

Recommended treatment duration: Begin with short sessions of 1–2 minutes and gradually increase treatment time. Most users benefit from sessions lasting 5–20 minutes.

Maximum treatment time: Do not exceed 20 minutes within a 6-hour period.

Recommended distance: Position yourself approximately 15–45 cm from the panel. Shorter distances provide higher light intensity, while greater distances allow treatment of a larger area.


Treatment area: Ensure the area being treated is fully exposed to the light and not covered by clothing or other materials.

Eye protection: Protective goggles are included. Avoid staring directly into the LEDs and always wear eye protection when using pulse mode.





Always consult a doctor before use in case of contraindications, medical conditions, or if you are unsure.

Note: The lamp does not emit heat! It is designed solely for LED therapy.


Product structure

What's included?


	Door hook Advanced Elite Multi		Hanging pulley system Advanced Elite Multi
	Hanging screws & cables Advanced Elite Multi		Remote controller Included with all models
	Protective eyeglasses Included with all models		Power cord Included with all models
	Wall bracket Ultra		

Accessories


Floor stand

Advanced
Elite
Multi


Floor stand

Ultra


Manual stand

Advanced
Elite

Electric stand

Advanced
Elite
Multi

Electric stand

Ultra


Installing the Door Hanger

STEP 1

Place the panel with the LED side facing toward you.


Attach the threaded mounting posts to the top of the panel by tightening them securely into place.

Then connect the suspension cables to the threaded posts.


STEP 2

Install the supplied door hook over a standard-sized door to create a secure hanging point for the panel.

STEP 3

Attach the braided suspension cables to the snap hook and carefully hang the panel from the door-mounted hook.


Installing the Door Hanger

STEP 4


To increase the height of the panel, support it with one hand while carefully pulling the adjustment cable with the other hand.

Important: Always keep the panel properly supported when adjusting its height using the pulley system.


STEP 5

To lower the panel, press the release lever on the height-adjustment pulley while securely supporting the panel with your free hand.


MULTI setup: connecting two Elite RLT panels

The PrimeFORCE Multi consists of two PrimeFORCE Elite panels that are connected together. Below, you can see how to assemble the panels easily and securely.

MULTI setup: connecting two Elite RLT panels

The panels can be connected using the Group function on the display.

You can find the instructions for this on page 23.

Touchscreen control panel


1 Home screen

This screen displays your current treatment settings, including the selected timer, wavelengths, brightness, and pulse frequency.

2 Presets

Select a preset program based on your wellness goals, such as skin health, fat loss, muscle recovery, joint support, sleep or pets.

To create a custom treatment, return to settings and manually adjust your preferred settings.

3 Settings

Access advanced controls to customize:

- Timer
- Mode
- Brightness
- Pulse
- Group Device
- Alarm clock
- Date & Time
- Other settings

4 Start / Pause

Press Start to begin your treatment session.

Press Pause to temporarily stop the session or end the treatment when desired.

Settings

Timer


Set treatment sessions from 1 to 30 minutes. You can adjust the timer in two ways:

Using the slider

Move the circular control around the timer ring to increase or decrease the session length.

Entering a time manually

Tap the center of the timer display, enter your preferred treatment duration, and confirm your selection.

Mode

RED


Activates the red light wavelengths: 630 nm, 650 nm, 660 nm, and 670 nm.

NIR (Near-Infrared)

Activates the near-infrared wavelengths: 810 nm, 830 nm, 850 nm, and 1060 nm.

RED + NIR

Activates all red and near-infrared wavelengths simultaneously: 630 nm, 650 nm, 660 nm, 670 nm, 810 nm, 830 nm, 850 nm, and 1060 nm.

Brightness


Adjust brightness from 0% to 100% to customize your treatment intensity. There are three ways to change the brightness level:

Using the “-” and “+” buttons

Press the “-” button to decrease brightness or the “+” button to increase brightness.

Using the slider

Move the brightness slider left or right until you reach your preferred intensity level.

Entering a value manually

Tap the brightness percentage, enter the desired value, and confirm your selection.

Pulse frequency

Adjust the pulse frequency anywhere from 0 Hz to 10,000 Hz. There are three ways to set the desired pulse rate:

Using the “-” and “+” buttons



Press the “-” button to lower the frequency or the “+” button to increase it.

Using the slider

Move the frequency slider until the desired pulse rate is displayed.

Entering a frequency manually

Tap the current pulse value, enter your preferred frequency, and confirm your selection.

Group

You can connect two panels using the “Group Device” function in the display settings.

Connect both panels to the power supply and turn them on. Then go to “Group Device” on both panels. Set panel A as “Primary” and panel B as “Secondary”.

Once the connection has been set up, you can control both panels using the panel that is set as “Primary”.

Alarm function

The built in alarm feature allows you to schedule a therapy session to start automatically. Go to "Alarm Clock" in the settings, tap "Edit" and add a new alarm using the plus icon.

The time you set is the end time of the treatment. For example, if you set 18:00 with a treatment duration of 10 minutes, the device will automatically start at 17:50 and stop at 18:00. This allows the treatment to run automatically before you wake up or during sleep, without needing to operate the device manually.

1 Alarm type

Select either a One-Time Alarm or a Daily Repeat Alarm, depending on your preference.

2 Configure Session Setting

Set the desired start time, session duration, wavelengths, brightness level, and pulse frequency for the alarm session

3 Manage Existing Alarms

Use the alarm menu to enable or disable alarms and edit their settings. Select an existing alarm to adjust the scheduled time or treatment parameters.

4 Add or Remove Alarms

The system supports multiple alarms. Use the "+" or "-" controls to add or remove alarms as needed, with support for up to 10 programmed alarms.

LIROMA

Hour Minute
01 04
02 05
03 06

Cancel OK

☐ only once ☒ Repeat

Alarm Time 08:00 >

Timer 10mins >

Mode RED/NIR >

Cancel OK

1

2

LIROMA

A 08:00 Repeat ☐ ○○○

01

A 08:00 Repeat ☐ ○○○

02

+ -

Run ↺

Group Device

Alarm Clock

Date&Time


Other Settings

3

4

Date & Time

The system clock can be adjusted to match your local date and time. Setting the correct date and time ensures that scheduled sessions and alarm functions operate as intended.

The image shows a digital screen with the LIROMA logo at the top. Below the logo is a time selection interface with columns for Hour, Minute, and Second. The current selection is 10:35:56. Below this are 'Cancel' and 'OK' buttons. Further down are four menu options: 'Group Device' (with a person icon), 'Alarm Clock' (with a clock icon), 'Date&Time' (with a clock icon), and 'Other Settings' (with a gear icon).

Hour	Minute	Second
9	34	55
10	35	56
8	36	57

Cancel OK

Group Device

Alarm Clock

Date&Time

Other Settings

Other settings

Work Alarm

When the work alarm is activated, the device will beep three times halfway through the treatment. For example, if the treatment duration is set to 10 minutes, the device will beep after 5 minutes.

Volume

Adjust the sound level of notifications or alarm signals

Screen Brightness

Adjust the screen brightness for comfortable use in different environments.

Language

Change the language of the control menu.

Factory Reset

Restore the panel to its original factory settings.

Red light therapy presets

Preset modes and recommended settings

Mode	630 nm	650 nm	660 nm	670 nm	810 nm	830 nm	850 nm	1060 nm	Timer	Pulse
Universal	100%	100%	100%	100%	100%	100%	100%	100%	15 min	None
Skin Care	80%	80%	0%	0%	80%	80%	0%	0%	10 min	None
Muscle	100%	100%	100%	100%	100%	100%	100%	100%	20 min	None
Recovery										
Joint Pain	100%	100%	100%	100%	100%	100%	100%	100%	20 min	None
Pets	50%	50%	50%	50%	100%	100%	100%	100%	20 min	None
Sleep	100%	100%	100%	100%	0%	0%	0%	0%	20 min	None
Fat Loss	100%	100%	100%	100%	100%	100%	100%	100%	20 min	None

How does it work?

● **Red light (630 nm, 650 nm, 660 nm, 670 nm)**

- Supports healthy skin function and appearance
- Stimulates collagen and elastin production
- Promotes circulation and nutrient delivery
- Supports cellular energy production (ATP)
- Helps improve skin texture and recovery

● **NIR (near-infrared light) (810, 830, 850, 1060 nm)**

- Penetrates deeper into muscles and joints
- Supports recovery and tissue regeneration
- Helps reduce inflammation and discomfort
- Promotes circulation in deeper tissues
- Supports cellular energy production at greater depths

Using RED + NIR together gives best results..

Understanding wavelengths

630 nm

The 630 nm wavelength falls within the visible red light spectrum and primarily works on the surface layers of the skin. This wavelength is commonly used for skin improvement, reducing redness, fine lines, acne, scars, and supporting collagen production. Because 630 nm is mainly active in the upper layers of the skin, it is particularly suitable for treatments focused on skin health, skin renewal, and achieving a fresher appearance.

660 nm

The 660 nm wavelength is one of the most commonly used wavelengths in red light therapy. It penetrates slightly deeper than 630 nm and supports both the skin and underlying tissues. 660 nm is often used to promote collagen production, wound healing, inflammation reduction, improved circulation, and recovery of superficial muscles and joints. This wavelength is also widely used for skin rejuvenation, fine lines, acne, and general cellular activity.

670 nm

The 670 nm wavelength is close to 660 nm but may penetrate slightly deeper. It is often associated with supporting the mitochondria, the energy-producing structures within cells. As a result, 670 nm may contribute to ATP production, or cellular energy. This makes the wavelength interesting for skin repair, muscle recovery, reducing inflammation, and supporting the body's natural healing processes.

810 nm

The 810 nm wavelength falls within the near-infrared spectrum and penetrates deeper than red light. It is commonly used for muscles, joints, tendons, nerve

tissue, and deeper healing processes. 810 nm is often considered beneficial for muscle soreness, joint discomfort, stiffness, and post-exercise recovery. This wavelength is also frequently associated with improved cellular energy production and support for the body's natural recovery capabilities.

830 nm

The 830 nm wavelength is a popular therapeutic wavelength within the near-infrared range. It penetrates deeply into tissue and is often used for muscle recovery, joint issues, tendons, connective tissue, and inflammation-prone areas. 830 nm may help improve circulation, support oxygen transport within cells, and reduce pain and stiffness. As a result, it is well suited for recovery, mobility, and support in cases of rheumatic conditions.

850 nm

The 850 nm wavelength is a deeply penetrating near-infrared wavelength. It is commonly used for deeper muscles, joints, and chronic inflammatory conditions. 850 nm supports mitochondrial activity, enabling cells to produce more energy. This may contribute to faster muscle recovery, reduced stiffness, improved circulation, and support for pain management. It is widely used in sports recovery, joint pain treatment, and general recovery therapies.

1060 nm


The 1060 nm wavelength penetrates even deeper than wavelengths in the 800 nm range. It is considered particularly interesting for deeper tissues, larger muscle groups, joints, and adipose (fat) tissue. In some applications, 1060 nm is used for deep heating effects, improved circulation, and supporting recovery processes in deeper tissues. This wavelength is also frequently mentioned in relation to body contouring and fat reduction, as it can penetrate deeper and interact with tissue differently than shorter wavelengths.

Summary

630, 660, and 670 nm are primarily suited for skin health, collagen production, redness reduction, wound healing, acne, and superficial inflammation.

810, 830, and 850 nm penetrate deeper and are mainly of interest for muscles, joints, tendons, circulation, pain relief, inflammation reduction, and recovery.

1060 nm penetrates the deepest and is particularly relevant for deeper tissues, larger muscle groups, joints, and potentially adipose (fat) tissue.

RLT Treatment Protocols

PrimeFORCE Series · 630–1060 nm · Multi-Wavelength Full Body Therapy

Muscle recovery, inflammation & sports performance

Red light (660 nm) and near-infrared (850 nm) reduce oxidative stress and inflammation in muscle tissue after intense exercise. The 850 nm wavelength penetrates up to 3–5 cm deep, allowing larger muscle groups to be effectively targeted.

Treat the stressed muscle groups before or immediately after training (within 2 hours). Hold the lamp at a distance of 15 cm and treat each zone for 3–5 minutes, corresponding to a target dose of 11–20 J/cm². Use the lamp 3 to 5 times per week and schedule 1 to 2 rest days. With consistent use, improvements in muscle recovery and reduction of muscle soreness (DOMS) are noticeable after 2 to 3 weeks.

- **Target dose: 10–20 J/cm²**
- **Time @ 15 cm: 3–5 min/zone**
- **Frequency: 3–5× per week**

Wound healing, skin improvement & anti-aging

660 nm red light stimulates fibroblast activity, promotes collagen synthesis, and accelerates cell regeneration. This makes it effective for both active wound healing and skin improvement/anti-aging applications. The difference lies in the dose: higher for wound healing, lower for anti-aging.

Wound healing

Hold the lamp perpendicular to the treatment area at a distance of 15 cm and treat for 90–120 seconds per zone (target dose 6–8 J/cm²). Treat daily or 5× per week for active wounds until healing has occurred.

- **Target dose: 6–8 J/cm²**
- **Time @ 15 cm: 90–120 sec/zone**
- **Frequency: daily or 5× per week**

Skin improvement & anti-aging

For collagen stimulation in healthy tissue: less is more. Treat for 30–60 seconds per zone (target dose 2–4 J/cm²). Consistency over weeks is more important than a higher dose per session. Visible results — such as improved skin texture and reduced fine lines — are typically noticeable after 4–8 weeks of regular use.

- **Target dose: 2–4 J/cm²**
- **Time @ 15 cm: 30–60 sec/zone**
- **Frequency: 5× per week**

Joint pain, inflammation & arthritis

Near-infrared light (850 nm) penetrates deep enough to reach the joint capsule and surrounding tissue. It inhibits pro-inflammatory cytokines (such as IL-1 β and TNF- α) and promotes local circulation, which can reduce pain and stiffness in both osteoarthritis and rheumatoid arthritis.

Direct the lamp at the joint from a distance of 15 cm. In the first two weeks, start with 4 minutes per zone (target dose ~15 J/cm²) and gradually increase to 8–13 minutes (~30–50 J/cm²) as the body adapts to the treatment. The full scientific therapeutic window extends up to 50 J/cm² (WALT guidelines; Bjordal et al., 2003) – higher doses within this range are particularly relevant for deeper joints such as hips or knees.

Treat 3 to 5 times per week; for chronic conditions, daily use during the first month is recommended. Results are typically noticeable after 3–6 weeks.

- **Target dose: 15–50 J/cm²**
- **Time @ 15 cm: 4–13 min/zone**
- **Frequency: 3–5× per week**

All protocols are based on published scientific research (including Leal-Junior et al. 2015, Avci et al. 2013, Brosseau et al. 2005, WALT guidelines). Consult a physician in case of serious or chronic conditions.

Warnings

Warnings and precautions

- Do not use this device on babies.
- Do not allow sleeping, disabled, or unconscious individuals to use this device.
- Do not use the device for more than 20 minutes per therapy session.
- Do not use the device on open wounds.
- Do not use the device in combination with liniments, ointments, or balms that contain heat-producing ingredients, as this may cause burns.
- Do not attempt to modify or repair the device yourself. If it is damaged or not functioning properly, please contact us as soon as possible at info@liroma.nl.

This device is electrical. To prevent electric shock and other hazards, please follow these instructions:

- Do not wash the product with water or other liquids.
- Do not place the device in water or in areas where it may be exposed to water or other liquids.
- Do not expose the device to running water.
- If the device comes into contact with or is submerged in water or any other liquid, immediately unplug it from the socket. If using the device near water, we recommend using a socket with a ground fault circuit interrupter (GFCI).
- Do not touch the device or its plugs with wet hands.
- Do not use or touch the device with wet or damp skin.
- Do not lift the device by the power cord.
- Do not use the device if the power cord is damaged.
- Do not use the device on or near heated surfaces.
- Do not use the device in areas where it may be exposed to flammable or combustible products or fumes, as this could cause a fire or explosion.
- Do not store the device at temperatures above 140°C (284°F).

Additional warnings

- Do not use the device if it is damaged, as continued use of a damaged device may cause electric shock or injury.
- Do not use extension cords unless they are designed to handle the total power of the device.
- Do not unplug the device by pulling on the cord. Always grasp the plug itself when removing it from the outlet.
- This device is not suitable for individuals with physical, neuro-sensory, or mental impairments, or those lacking experience and/or knowledge.
- Do not use the lamp in combination with implanted electronic devices, such as artificial hearts or artificial lungs, without first consulting your doctor.
- Do not use the lamp on pregnant women, individuals with severe heart conditions, critically ill patients, or emergency cases without consulting your doctor.
- Discontinue treatment immediately if you experience hypersensitivity or allergic skin reactions due to medication, cosmetics, or food.
- Under no circumstances should you open the device or attempt to repair it yourself, as this may cause injury and void your warranty. If the device is damaged or not functioning properly, contact the retailer or info@liroma.nl as soon as possible.

Do not use this product without first consulting your doctor if you suffer from:

- Severe illnesses such as heart disease, acute inflammation, thrombosis, blood clotting disorders, adrenal suppression, systemic lupus erythematosus, or malignant conditions.
- Hypersensitivity to infrared light or have edema.
- A condition related to nerve demyelination, such as carpal tunnel syndrome, multiple sclerosis, or ulnar neuropathy.
- Pregnancy.
- Certain medications or plant extracts that cause photosensitivity are not compatible with light therapy. (Examples: antibiotics, antidepressants, Roaccutane, MTX, antifungal medication, and St. John's Wort).

Important notes

- Red light is visible, while NIR infrared light is invisible to the naked eye. However, with NIR infrared light, you may perceive a slight glow or dimming, which is completely normal.

Care and cleaning

Follow the instructions below to clean and maintain the device properly and safely.

- Turn off and unplug
- Let device cool
- Clean with dry cloth
- Do not use liquids
- Do not open device

Liroma has been a specialist in red light therapy for years. We offer lamps that provide full-body exposure, as well as smaller devices that allow for targeted treatment on specific body parts. Additionally, we have red light devices specifically designed for the back and face. Visit our website to explore all our products.

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
Post a review with
a photo and receive
€10 cashback.

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
The screen displays the LIROMA logo at the top. Below it, there are two columns for 'Hour' and 'Minute'. The 'Hour' column shows '01', '02' (highlighted), and '03'. The 'Minute' column shows '04', '05' (highlighted), and '06'. Below these columns are 'Cancel' and 'OK' buttons. Further down, there are two radio buttons: 'only once' and 'Repeat' (selected). Below these are three rows of settings: 'Alarm Time' set to '08:00 >', 'Timer' set to '10mins >', and 'Mode' set to 'RED/NIR >'. At the bottom are 'Cancel' and 'OK' buttons.

The screen displays the LIROMA logo at the top. Below it, there are two rows of alarm settings. Each row shows a label 'A' followed by a time '08:00' and a 'Repeat' toggle switch. The first row is labeled '01' and the second '02'. Below these are two large circular buttons: a '+' button and a '-' button. Below these are two buttons: 'Run' and a button with a circular arrow icon. At the bottom, there are four buttons: 'Group Device', 'Alarm Clock', 'Date&Time', and 'Other Settings'.

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